

**Arikaree Environmental Laboratory**

**Revision 1**

**Date: Oct. 10, 2018**

**NO<sub>3</sub>+NO<sub>2</sub>(N) by FIA SOP**

**TITLE: Lachat QuikChem Method 10-107-04-1-A**

**Lachat QuikChem Method 10-117-07-1-B**

**ANALYTES: NO<sub>3</sub>+NO<sub>2</sub>(N)**

---

**1) Applicable Matrices**

- a) This method is applicable to ground, surface and waste waters.
- b) The method is applicable in the range from 0.01 to 10.0 mg/L for NO<sub>3</sub>+NO<sub>2</sub>(N) analysis.

**2) Scope and Application**

- a) Nitrate is quantitatively reduced to nitrite by passage of the sample through a copperized cadmium column. The nitrite produced is then diazotized with sulfanilamide followed by coupling with N-(1-naphthyl)ethylenediamine dihydrochloride (NED). The resulting complex has a magenta color and its absorbance at 520nm is quantified.

**3) Interferences**

**NO<sub>3</sub>+NO<sub>2</sub>(N)**

- a) Residual chlorine can interfere by oxidizing the cadmium reduction column.
- b) High concentrations of iron, copper and other metals can create biased low results. EDTA is added to the ammonium chloride buffer to complex with the metals and reduces this interference.
- c) Samples containing high amounts of oil and grease can coat the surface of the reduction column.
- d) Turbid samples may interfere with light absorption. Filtering samples through a 0.45µm membrane filter before analyzing will remove turbidity.

**4) Equipment and Supplies**

- a) Analytical balance capable with 0.0001g sensitivity.
- b) Top loading balance with 0.01g sensitivity.
- c) Pipettes capable of precisely dispensing 1.0-5.0mL and 0.1-1.0mL
- d) Pipette tips appropriate for each pipette.
- e) Class A volumetric flasks and pipettes for reagent and standard preparation.

- f) Lachat 8000 flow-injected autoanalyzer with  $\text{NO}_3+\text{NO}_2(\text{N})$  manifold, 10.00mm flow cell, and interference filter at 520 nm.
- g) Water purifier for preparation of ASTM Type I water (Ultrapure).
- h) Various Class A volumetric pipettes for dilution of calibration standards.

## 5) Reagents and Standards

### $\text{NO}_3+\text{NO}_2(\text{N})$

- a) 10N NaOH: In a volumetric flask filled with ~400mL of ultrapure water, dissolve 400g of NaOH pellets. Stir until pellets are dissolved, adding water if necessary. Allow solution to cool (place in water bath) and dilute to the mark. Store in a 1L Nalgene and store in the base cabinet.
- b) Ammonium Chloride buffer: In a class A 1L volumetric flask dissolve 85g ammonium chloride, 1.0g ethylenediamine tetraacetic acid disodium salt dihydrate and 11.0mL 10N NaOH. Dilute to mark.
- c) Sulfanilamide Color reagent: In a class A 1L volumetric flask dilute 100mL 85% phosphoric acid, dissolve 1.0g N-(1-Naphyl)ethylenediamine Dihydrochloride (NED) and 40.0g sulfanilamide. Dilute to the mark.

## 6) Sample Collection, Preservation, Shipment, and Storage

- a) Collect and filter samples in clean plastic bottles. Samples may be analyzed unpreserved, but sulfuric acid may be added to preserve the samples at <pH 2, or the samples may be stored frozen for preservation.
- b) Acid preserved samples may be stored up to 28 days. Frozen samples more be stored indefinitely.

## 7) Quality Control

- a) Calibrate before each analytical run. Analyze Instrument Performance Check (IPC)/ from a second source standard after calibration. Calibration r-value is required to be > 0.995, if not re-calibrate.
- b) Run a blank (ICVB) at the beginning of each run and then once every 20 samples. ICVB must be <LOD. If not, re-analyze, if still out of control qualify data. Subsequent blanks (CCVB) must be <LOD, or less than 10% of the measured concentration in the samples from the adjacent set.
- c) A 5.0 mg/L Nitrite standard must be analyzed after the blank along with a 5.0 mg/L Nitrate standard. The nitrate standard can be an ICVS. Calculate the reduction column efficiency using the ratio of the response factors of nitrate/nitrite. The efficiency must be  $\geq 75\%$ . If below this, replace/regenerate the column.
- d) Analyze a Continuing Calibration Verification Standard (CCVS) with every 20 samples.

- e) A duplicate and matrix spike must be analyzed every 20 samples (5%). Prepare a 50mg/L matrix spike solution as follows: Pipet 5.0mL of 1,000mg/L NO<sub>3</sub><sup>-</sup> stock standard into a class A 100mL volumetric flask. Dilute to the mark, replace the stopper and invert to mix. Store in a 100mL Nalgene bottle.

## 8) Calibration

- a) Prepare high range standards as listed below. Use only class A volumetric pipettes and flasks. Dilute using ultrapure water.

Standard	NO <sub>3</sub> Concentration (mg/L)	mL NO <sub>3</sub> Standard	Total volume
1	10.0	5.0	500mL
2	5.0	5.0	1,000mL
3	1.0	50.0 of Standard 1	500mL
4	0.5	50.0 of Standard 2	500mL
5	0.1	5.0 of Standard 1	500mL
6	0.05	5.0 of Standard 2	500mL

- b) Prepare a 5.0 mg/L Nitrite standard. Fill a 1,000mL class A volumetric flask with ~500mL of ultrapure water. Pipet 5.0mL, using a 5.0mL class A pipet, of 1,000ppm Nitrite stock standard into the 1,000mL flask. Dilute to the mark on the flask with ultrapure water. Invert to mix

## 9) Procedure for Lachat 8500 operation

- a) Turn on the Lachat autoanalyzer by moving the switch on the power strip behind the autosampler to the on position. The instrument will go through its normal startup routine (autosampler probe moves, switching valves on the instrument turn, etc.).
- b) Place reagent tubing into water and start the peristaltic pump. Clamp each pump tube clamp onto the rollers. Push the tensioner cam lever back one click.
- c) Hit the space bar on the computer keyboard to “wake-up” the computer. Double-click on the Omnion icon on the desktop.
- d) Click the open button on the top toolbar. Find the correct template for NO<sub>3</sub>+NO<sub>2</sub>(N) analysis and double-click to open.
- e) Make sure the concentrations for the IPC/QCS sample are correct by clicking on the sample in the worksheet. If the values are incorrect, change them by highlighting all of the samples in the DQM set (most likely DQM1) and clicking “clear DQM set” from the menu. Change the concentration in the “Run Properties” window in the upper right of the screen. Highlight the DQM samples again and right click. Select define DQM set. Set the frequency to once.

- f) Start inputting sample names/numbers into the worksheet. To add more spaces for samples, highlight a row and right click. Select "Insert many" and input how many rows to add. Remember to hit "enter" after each sample is input. Use the benchsheet to enter the sample identifiers. To enter multiple samples quickly, select multiple rows and right click. Go to "Columns" and select "Auto SampleID." Enter the part of the sample identifier that will not change in the "Fixed Part of Sample ID" then the number or letter to start. Click "Accept." It is helpful if the number of rows highlighted matches the number of samples to be input automatically.
- g) Once finished, highlight the entire worksheet and right click. Go to "Columns" and select "Auto-number cups." Make sure the number of cups in the worksheet matches what is on the benchsheet. If the numbers do not match, check the duplicates and matrix spikes to see if they match. Once you find a matching duplicate and spike set, you will know where the discrepancy is located.
- h) Remove the reagent lines from the water and begin placing them in their corresponding reagent bottles. Allow reagents to flow through manifold for ~10 minutes.
- i) Click on the "Configuration" menu at the top of the Omnion window. Go to "Autosamplers" and click the button "Intialize Autosampler." This will put the sampler probe into the wash well to begin drawing water through the sample line.
- j) After reagent has been pumping for ~10 min. click the preview button at the top of the Omnion window. This will let you see the baseline before starting the run. The baseline is flat and the signal is where it is expected, click "Stop" next to the "Preview" button and then click "Start."
- k) During the beginning of the run, watch to make sure there is good Gaussian peak shape and the peak expection window is integrating the entire peak correctly.
- l) After the calibration standards are integrated, check the calibration curve by clicking on the bottom button on the left side of the channel window. Make sure calibration coefficient is  $\geq 0.995$  and %RSD for each standard is under 10%. If either of these limits are exceeded begin by excluding peaks that may skew the curve. If this does not correct the problem, recalibrate.
- m) Next, be sure the initial calibration check standards are within 10% of the expected value. If not, reanalyze the standards. If this does not correct the issue, recalibrate.
- n) After everything in the beginning is in control, the rest of the run can be monitored for duplicate RPD's and matrix spike recoveries. If these are out of control, try reanalyzing the samples. Also, a continuing calibration verification blank (CCVB) and continuing calibration verification standard will be analyzed every 20 samples. Make sure these are within 10% of known value.
- o) If everything in the run goes to plan, the instrument will stop after all of the samples in the worksheet have been analyzed. The reagent lines can be removed and placed into ultrapure water for ~10min. After the 10min., remove tubes from water and hang them on the hook above the pump and allow air to be pumped through the manifold. Pump air until all water has been evacuated from the manifold lines.

- p) Release the pump tubes on the pump and stop the pump. Turn of the Lachat by flipping the switch on the power strip behind the autosampler.
- q) Samples in the test tubes can be disposed by dumping the sample into the waste drum. The tubes will be cleaned and used again.

#### **10) Corrective Action for Out-of-Control Data**

- a) Any CCV blank that exceeds the LOD, the analyst must inspect the concentration of the previous sample. If the blank is greater than 10% of the previous sample, reanalyze the blank immediately following the CCVS. No samples can be run until the blank meets requirements.
- b) CCV and LCS standards must fall within 10% of true value (90%-110% Recovery). If not re-mix and re-analyze, if still out of range re-calibrate. Fill out corrective action report.

#### **11) Contingencies for Handling Out-of-Control Data**

- a) Samples that fail the CCVB or ICVS will have to be qualified back to the last sample that the quality control met the above conditions.

#### **12) Waste Management**

- a) All waste must be collected (except for waste from probe wash well) and disposed of in the correct hazardous waste drum.

#### **13) References**

- a) Standard Methods for the Examination of Water and Wastewater, 22nd ed. Clesceri, L.S.; Greenberg, A.
- b) Lachat QuikChem Method 10-107-04-1-A. Determination of Nitrate/Nitrite in Surface and Wastewaters by flow injection analysis. K. Wendt, revision date 25 March 2003.